




Improve Your Granulator's Overall Performance and Save Money Too!

Stationary Bed Knives and Rotating Fly Knives

- Stationary knife blade can be sharpened individually
- Rotating knife blades should be sent in sets for matched grinding and balancing
- Specify machine make and type
- Specify any minimum resharpening dimensions
- Quick turnaround



Fly Knife



Bed Knife

Conair/Wortex Fly Knives

- Diameter of the cutting circle remains unchanged
- Special diamond grinding tools sharpen the underside of insert only
- Includes sharpening of all 3 inserts



BEFORE

Segmented Fly Knives (3-tooth)
\$6.75/Tooth



AFTER

For New Replacement Granulator Knife Blades, call us at **800-627-1033** or visit www.servicesforplastics.com

EXAMPLES:

